



Multitron Standard



INFORS 
We bring life to your laboratory

This is the Multitron Standard

This incubation shaker comes standard with everything you need for microbial applications and insect cell cultures. All you need to determine is the number of units, the shaking throw and if a cooling is needed—and then you can get right to work.

Applications

- Suitable for microorganisms and insect cells, e.g., mRNA, protein expression or media development.

Easy to use

- Automatic restart after closing the door and automatically locking, extendable trays make daily laboratory work easier.

Comfortable working height

- The 31 cm tall base provides comfortable access when working with one and two units. When three units are stacked on top of each other, **the maximum working height is just 1.3 m.**

Reproducible results

- High-quality components and the latest technology ensure reliable, long-term operation. Uniform temperature distribution across the entire tray guarantees homogeneous conditions and thus reproducible and comparable results.

Versatile configuration options

- Available as a single unit or as two or three units stacked on top of each other; you may expand your system at any time.

Large capacity

- Despite their small footprint of just under one square meter, three stacked Multitron Standard units allow you to incubate over 60 liters of culture liquid. The interior can accommodate Erlenmeyer flasks up to 5 L in size.

	Base	Cooling	Shaking throw
Single unit	31 cm	In base	25/50 mm
Single unit	13 cm or 31 cm	Top cooling	25/50 mm
Two stacked units	13 cm or 31 cm	Top cooling	25/50 mm
Three stacked units	13 cm	Top cooling	25/50 mm



Multitron Standard,
three stacked units
with top cooling,
working height: 1.30 m

Up to 6 x 5 L
Erlenmeyer flasks

Compact footprint
of just under one
square meter

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

