



Whitley Workstations Single Plate Entry Systems



**Introduce individual samples
without compromising internal
atmospheric conditions**

**Versions for the majority of our
workstations can be retro fitted**

**Perfect for on-call and emergency
use**

The Whitley Workstation range includes a variety of options and accessories enabling you to tailor a system to your needs and your laboratory working practices.

A very useful option is a single plate entry system, which allows individual Petri dishes to be introduced into the chamber without compromising internal atmospheric conditions. Multiwell plates, cell culture flasks and similar sized objects can be introduced using the same system.

The single plate entry system can be extremely handy for on-call and emergency use where individual samples can be introduced quickly without the need to use the airlock or ports.

More convenient working methods and lower gas usage, what more could you want?

	9cm*				15cm	
	A03125	A03225	A03325	A07208	A06103	
DG250	●	●	●	-	-	DG250
A20	-	-	-	-	●	A20
A25	-	-	-	-	●	A25
A35	-	-	-	●	●	A35
A45	-	-	-	●	●	A45
A55	-	-	-	●	-	A55
A85	●	●	●	-	-	A85
A95	●	●	●	-	-	A95
A135H	-	-	-	●	●	A135H
A135GMP	-	-	-	-	-	A135GMP
A155H	-	-	-	●	-	A155H
A155GMP	-	-	-	-	-	A155GMP
M35	-	-	-	●	●	M35
M45	-	-	-	●	●	M45
M85	●	●	●	-	-	M85
M95	●	●	●	-	-	M95
M135H	-	-	-	●	●	M135H
M155H	-	-	-	●	-	M155H
H25	-	-	-	-	●	H25
H35	-	-	-	●	●	H35
H45	-	-	-	●	●	H45
H85	●	●	●	-	-	H85
H95	●	●	●	-	-	H95
H135H	-	-	-	●	●	H135H
H135GMP	-	-	-	-	-	H135GMP
H155H	-	-	-	●	-	H155H
H155GMP	-	-	-	-	-	H155GMP
i2	-	-	-	-	-	i2
	w: 107mm h: 26mm Front loading Can be retro fitted	w: 107mm h: 26mm Front loading Longer tray for larger items Can be retro fitted	w: 98mm h: 20mm Front loading With folding tray for an unobstructed work area Can be retro fitted	w: 107mm h: 26mm Front loading Can be retro fitted on units with a removable front	w: 152mm h: 52mm Located on the left hand panel of the unit Can be retro fitted	

Key: - Not available / ● Option available

Single Plate Entry Systems are designed for occasional use and not as a direct alternative to airlocks or walk-in rooms that allow personnel to pass through to a clean room. Detailed information for installation

Single Plate Entry Systems/004-05/March 2021

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

